

Dimethyl silicone oil (viscosity 100 cSt (25 °C))

Chemical Properties

CAS No. :

Formula: $[\text{Si}(\text{CH}_3)_2\text{O}]_n$

Molecular Weight:

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dimethyl silicone oil, viscosity 100 cSt (25 °C) (Poly(dimethylsiloxane) (PDMS), viscosity 100 cSt (25 °C); Dimethylsilicone fluid, viscosity 100 cSt (25 °C)) is a biochemistry reagent.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481