

Polyvinylpyrrolidone K120

Chemical Properties

CAS No. :

Formula: (C₆H₉NO)_n

Molecular Weight:

Storage: Powder: -20°C for 3 years
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Polyvinylpyrrolidone K120 is extensively utilized in both human and veterinary medicine as an effective wound healing promoter and disinfectant when combined with iodine and other compounds.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481