

mPEG-amine (MW 20000)

## Chemical Properties

CAS No. :

Formula: (C<sub>2</sub>H<sub>4</sub>O)<sub>n</sub>C<sub>3</sub>H<sub>9</sub>NO

Molecular Weight:

Storage: Keep away from direct sunlight, Store under nitrogen  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	m-PEG-NH <sub>2</sub> (MW 20000) is a PROTAC linker that falls under the PEG category. It can be utilized in the synthesis of PROTAC molecules.
Targets(IC <sub>50</sub> )	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481