

mPEG-amine (MW 10000)

Chemical Properties

CAS No. :

Formula: (C₂H₄O)_nC₃H₉NO

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	mPEG-amine (mPEG-NH ₂) (MW 10000) is a PROTAC linker from the PEG category, employed in synthesizing PROTAC molecules.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481