

Poly(acrylic acid) (MW 450000)

## Chemical Properties

CAS No. :

Formula: (C<sub>3</sub>H<sub>4</sub>O<sub>2</sub>)<sub>x</sub>

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Poly(acrylic acid) (MW 450000) is a 450000 molecular weight anionic polymer. It serves as a corrosion inhibitor and surface stabilizer.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481