

PLGA (75:25)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PLGA (75:25) is a low-toxicity carrier with biocompatibility and biodegradability, designed for the controlled release of active molecules and achieving sustained release within biological systems. Composed of 75% polylactic acid (PLA) and 25% polyglycolic acid (PGA), PLGA (75:25) has been extensively utilized in the research of delivering active molecules, proteins, and various macromolecules such as DNA, RNA, and peptides.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481