

Poly-L-lysine hydrobromide (MW 4000-15000)

## Chemical Properties

CAS No. :

Formula: C<sub>8</sub>H<sub>19</sub>BrN<sub>2</sub>O

Molecular Weight:

Keep away from moisture

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

Poly-L-lysine (hydrobromide) (MW 4000-15000) enhances cell adhesion to solid substrates. It is utilized as a biomaterial or organic compound in life sciences research.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481