

Poly-L-lysine hydrobromide (MW 15000-30000)

Chemical Properties

CAS No. :

Formula: C₈H₁₉BrN₂O

Molecular Weight:

Keep away from moisture

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Poly-L-lysine (hydrobromide) (MW 15000-30000) enhances the adhesion of cells to solid-phase substrates. It can be utilized as a biomaterial or an organic compound in life science research.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481