

Poly-L-lysine solution (0.01%, sterile-filtered, MW 150000-300000)

## Chemical Properties

CAS No. :

Formula: (C<sub>6</sub>H<sub>14</sub>N<sub>2</sub>O<sub>2</sub>)<sub>n</sub>.nBrH

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

Poly-L-lysine solution (0.01%, sterile-filtered, MW 150000-300000) (PLL solution (0.01%, sterile-filtered, hydrobromide, MW 150000-300000, for cell culture)) is a biochemistry reagent for cell culture.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481