


FT385

Chemical Properties

CAS No. : 2414580-72-4
 Formula: C₂₂H₁₉N₅O₂
 Molecular Weight: 385.42
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	FT385 is a novel USP30 inhibitor, recapitulating genetic loss of USP30 and setting the trigger for PINK1-PARKIN amplification of mitochondrial ubiquitylation.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5946 mL	12.9729 mL	25.9457 mL
5 mM	0.5189 mL	2.5946 mL	5.1891 mL
10 mM	0.2595 mL	1.2973 mL	2.5946 mL
50 mM	0.0519 mL	0.2595 mL	0.5189 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481