



Anti-BAX Polyclonal Antibody 4

Product Details

Ig Type: IgG
Reactivity: Human, Mouse, Rat, Rabbit (predicted: Dog, Pig, Cow, Sheep)
Molecular Weight: Theoretical: 21 kDa.
Purification: Protein A purified

Applications

1. Blank control (blue line): HeLa (fixed with 80% methanol (5 min at -20°C) and then permeabilized with 0.1% PBS-Tween for 20 min at room temperature).
Primary Antibody (green line): Rabbit Anti-Bax antibody (TMAB-00191), Dilution: 1 µg/10⁶ cells;
Isotype Control Antibody (orange line): Rabbit IgG.
Secondary Antibody (white blue line): Goat anti-rabbit IgG-FITC, Dilution: 1 µg/test.
Verified Activity:
2. Sample:
HeLa (Human) Cell Lysate at 30 µg
Primary: Anti-Bax (TMAB-00191) at 1/1000 dilution
Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution
Predicted band size: 21 kDa
Observed band size: 23 kDa



Application: FCM, ICC/IF, IHC-Fr, IHC-P, WB
Recommended: WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; ICC/IF: 1:100; IF: 1:100-500; FCM: 1 µg/test

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: KLH conjugated synthetic peptide: human Bax

Antigen Species: Human

Gene ID: 581

Uniprot ID: Q07812

Synonyms: BCL2L4;Bcl-2-like protein 4;Bcl2-L-4;Apoptosis regulator BAX

Biology Area: Bcl 2 family, Metabolism, Bcl 2 family, Associated Proteins, Bcl2 Family, Apoptosis marker and proteins ELISA kits, Apoptosis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481