


## Beta-Actin (Fruit Fly) Loading Control Polyclonal Antibody

## Product Details

Ig Type: IgG  
Reactivity: Fruit Fly  
Molecular Weight: Theoretical: 42 kDa. Actual: 42 kDa.  
Purification: Protein A purified

## Applications

Sample:  
Lane 1: Fruit Fly lysates  
Lane 2: Fruit Fly lysates  
Verified Activity: Primary: Anti-Actin (Fruit Fly) (TMAB-00213) at 1/1000 dilution  
Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution  
Predicted band size: 42 kDa  
Observed band size: 42 kDa



Application: WB  
Recommended WB: 1:1000-5000

## Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

## Antigen Details

Immunogen: Recombinant Protein: Fruit Fly Actin.  
Synonyms: Dro-actin;Act5C; $\beta$ -actin;ACT1\_DROME;Fruit Fly Actin;Actin 5C  
Biology Area: Beta Actin,Actin,Actin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481