

Anti-Beta-Amyloid 1-40 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	(predicted:Human,Mouse,Rat,Chicken,Dog,Pig,Cow,Rabbit)
Molecular Weight:	Theoretical: 4.4 kDa.
Purification:	Protein A purified

Applications

Application:	ELISA,ICC/IF,IF,IHC-Fr,IHC-P
Recommended	IHC-P: 1:100-500; IHC-Fr: 1:100-500; ICC/IF: 1:100-500; IF: 1:100-500; ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human beta-Amyloid
Antigen Species:	Human
Gene ID:	351
Uniprot ID:	P05067
Synonyms:	APP695;Amyloid Precursor;Amyloid Precursor Protein 695
Biology Area:	Regulation,RNA Processing,DNA / RNA binding

Research Background

The cerebral and vascular plaques associated with Alzheimer's disease are mainly composed of Amyloid beta peptides. beta Amyloid is derived from cleavage of the Amyloid precursor protein and varies in length from 39 to 43 amino acids. beta Amyloid [1-40], beta Amyloid [1-42], and beta Amyloid [1-43] peptides result from cleavage of Amyloid precursor protein after residues 40, 42, and 43, respectively. The cleavage takes place by gamma-secretase during the last Amyloid precursor protein processing step. beta Amyloid [1-40], beta Amyloid [1-42], and beta Amyloid [1-43] peptides are major constituents of the plaques and tangles that occur in Alzheimer's disease. beta Amyloid antibodies and peptides have been developed as tools for elucidating the biology of Alzheimer's disease.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481