

Anti-Beta-lactoglobulin Polyclonal Antibody 2

Product Details

Ig Type:	IgG
Reactivity:	(predicted:Cow,Sheep)
Molecular Weight:	Theoretical: 18 kDa.
Purification:	Protein A purified

Applications

Application:	ELISA,IF,IHC-Fr,IHC-P,WB
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500; ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: cow Beta-lactoglobulin
Antigen Species:	Cow
Gene ID:	280838
Uniprot ID:	Q32P89
Synonyms:	Bos d 5;Lactoglobulin beta;beta-LG;progestagen-associated endometrial protein;BLG;LGB;beta lactoglobulin B;beta-lactoglobulin precursor;LACB_BOVIN;β Lactoglobulin
Biology Area:	Lipid metabolism

Research Background

bov beta-lactoglobulin is the content of the milk in the category of a protein myosin. Protein contained in whole milk (3.7 percent) to 80% of the casein in the removal of the casein whey protein, 55% of the beta-lactoglobulin. Access to the crystallization of the isoelectric point = 5.2, the pH in this dual capacity under the main structure. The people in latex-free Beta-lactoglobulin. The beta-lactoglobulin and bov beta-lactoglobulinA-based.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481