

Anti-BVDV E2 Polyclonal Antibody

Product Details

Ig Type: IgG
Reactivity: (predicted:BVDV)
Molecular Weight: Theoretical: 45 kDa.
Purification: Protein A purified

Applications

Application: ELISA
Recommended ELISA: 1:5000-10000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinded BVDV Envelope glycoprotein E2
Synonyms: PestiV1gp1 polyprotein [Bovine viral diarrhea virus 1];POLG_BVDVN;Envelope glycoprotein E2;
BVDV Envelope glycoprotein E2;polyprotein [Bovine viral diarrhea virus 1]

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481