

## Anti-Flagellin B(flxB) Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	(predicted:Vibrio Parahemolyticus)
Molecular Weight:	Theoretical: 40 kDa.
Purification:	Protein A purified

### Applications

Application:	ELISA,ICC/IF,IF,IHC-Fr,IHC-P,WB
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; ICC/IF: 1:100-500; IF: 1:100-500; ELISA: 1:5000-10000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	KLH conjugated synthetic peptide: Vibrio Parahemolyticus Polar flagellin B/D
Synonyms:	flxB;Flagellin B

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481