


Anti-FMDV VP1 (type O) Antibody (7B798)

Product Details

Ig Type:	IgG1
Reactivity:	FMDV
Molecular Weight:	Actual: 46 kDa.
Clone:	7B798
Purification:	Protein G purified

Applications

Sample:	Lane 1: Recombinant FMDV VP1 (type O) protein, Sumo & His
Verified Activity:	Primary: Anti-FMDV VP1 protein (type O) (TMAB-00692) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution Observed band size: 46 kDa



Application:	ELISA, WB
Recommended	ELISA=1:5000-10000; WB=1:500-2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombined FMDV VP1 protein (type O)
Synonyms:	Foot-and-mouth disease virus-type O; FMDV VP1; FMDV Polyprotein(VP1 protein); VP1 protein; foot and mouth disease virus

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481