


Anti-FMDV VP1 (type O) Polyclonal Antibody

Product Details

Ig Type: IgG
Reactivity: FMDV
Molecular Weight: Actual: 46 kDa.
Purification: Protein A purified

Applications

Sample:
Lane 1: Recombinant FMDV VP1 (type O) protein, Sumo & His
Verified Activity: Primary: Anti-FMDV VP1 protein (type O) (TMAB-00693) at 1/1000 dilution
Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution
Observed band size: 46 kDa



Application: ELISA, WB
Recommended: ELISA=1:5000-10000; WB=1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombined FMDV VP1 protein (type O)
Synonyms: VP1 protein; foot and mouth disease virus; Foot-and-mouth disease virus-type O; FMDV Polyprotein(VP1 protein); FMDV VP1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481