

Anti-HBeAg Antibody (5E79)

Product Details

Ig Type:	IgG
Reactivity:	HBV
Molecular Weight:	Theoretical: 18 kDa. Actual: 25 kDa.
Clone:	5E79
Purification:	Protein G purified

Applications

Sample:

Lane 1: Recombinant HBcAg protein, Tag free

Lane 2: Recombinant hepatitis B virus e virus like particle protein, His


Verified Activity:

Primary: Anti-HBeAg (TMAB-00832) at 1/1000 dilution

Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution

Predicted band size: 18 kDa

Observed band size: 25 kDa



Application: ELISA, WB

Recommended ELISA=1:500-5000; WB=1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen:	HBeAg full length protein
Synonyms:	Precore protein; HBEAG_HBVD3;p25;C;External core antigen
Biology Area:	Hepatitis B

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481