

Anti-HBeAg Antibody (7X515)

Product Details

Ig Type:	IgG
Reactivity:	(predicted:HBeAg)
Molecular Weight:	Theoretical: 18 kDa.
Clone:	7X515
Purification:	Protein G purified

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	HBeAg full length protein
Synonyms:	External core antigen;Precore protein;C;HBEAG_HBVD3;p25
Biology Area:	Hepatitis B

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481