

Anti-HBsAg Antibody (8A718)

Product Details

Ig Type:	IgM
Reactivity:	(predicted:HBsAg)
Molecular Weight:	Theoretical: 44 kDa.
Clone:	8A718
Purification:	Affinity purified

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: HBsAg protein
Synonyms:	LHB;Large S protein;L glycoprotein;HBSAG_HBVC5;Large surface protein;Major surface antigen;S;L-HBsAg
Biology Area:	Hepatitis B

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481