

Anti-HPV16 E6 Antibody (3K405)

Product Details

Ig Type:	IgG
Reactivity:	HPV16
Molecular Weight:	Theoretical: 11 kDa.
Clone:	3K405
Purification:	Protein A purified

Applications

Application:	ELISA,IF,IHC-Fr,IHC-P
Recommended	IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500; ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: HPV16 E6 protein
Synonyms:	VE6_HP16;HPV16 E6;Human Papilloma Virus;Human papillomavirus type 16;Human papillomavirus type 16 E6;E6 protein;E6;HPV16-E6;Protein E6

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481