

Anti-HPV16 E7 Polyclonal Antibody 4

Product Details

Ig Type: IgG
Reactivity: (predicted:HPV16)
Molecular Weight: Theoretical: 11 kDa.
Purification: Protein A purified

Applications

Application: FCM,IF,IHC-Fr,IHC-P
Recommended IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500; FCM: 1ug/Test

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: KLH conjugated synthetic peptide: HPV16-E7 protein
E7 protein(HPV16);HPV16 E7 protein;Human Papilloma Virus;VE7_HP16;HPV16 E7;Protein E7;
Synonyms: Human papillomavirus type 16 E7;HPV 16 E7 protein;E7;Human Papillomavirus 16(E7);HPV16-E7

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481