

## Anti-KLH/Keyhole Limpet Hemocyanin Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	(predicted:Keyhole limpet hemocyanin (KLH))
Molecular Weight:	Theoretical: 800-900 kDa.
Purification:	Protein A purified

### Applications

Application:	ELISA
Recommended	ELISA: 1:5000-10000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Full length protein: KLH.
Synonyms:	KLH 1;Keyhole limpet hemocyanin(KLH);Hemocyanin 1;Keyhole limpet hemocyanin A;KLH-A;hemocyanin A chain;KLH1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481