

Anti-LTBP4 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	(predicted:Human,Mouse,Rat,Dog,Pig,Cow,Sheep)
Molecular Weight:	Theoretical: 170 kDa.
Purification:	Protein A purified

Applications

Application:	WB,IHC-P,IHC-Fr,IF,ELISA
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500; ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human LTBP4
Antigen Species:	Human
Gene ID:	8425
Uniprot ID:	Q8N2S1
Synonyms:	Ltbp4;LTBP-4;Latent-transforming growth factor beta-binding protein 4
Biology Area:	TGF,TGF,Other ECM Proteins,TGF

Research Background

May be involved in the assembly, secretion and targeting of TGFB1 to sites at which it is stored and/or activated. May play critical roles in controlling and directing the activity of TGFB1. May have a structural role in the extra cellular matrix (ECM).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481