

Anti-Meteorin-like/METRNL Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	(predicted:Human,Mouse,Rat,Dog,Cow,Horse,Sheep)
Molecular Weight:	Theoretical: 30 kDa.
Purification:	Protein A purified

Applications

Application:	ELISA,ICC/IF,IHC-Fr,IHC-P
Recommended	ELISA=1:5000-10000; ICC/IF=1:100-500; IHC-Fr=1:100-500; IHC-P=1:100-500

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human METRNL
Antigen Species:	Human
Gene ID:	284207
Uniprot ID:	Q641Q3
Synonyms:	METRNL;Subfatin;Meteorin-like protein
Biology Area:	Other Antibodies

Research Background

METRNL is a 311 amino acid secreted protein belonging to the Meteorin family and may have similar roles to that of the Meteorin protein. METRNL is encoded by a gene located on human chromosome 17, which comprises over 2.5% of the human genome and encodes over 1,200 genes, some of which are involved in tumor suppression and in the pathogenesis of Li-Fraumeni syndrome, early onset breast cancer and a predisposition to cancers of the ovary, colon, prostate gland and fallopian tubes.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481