

Anti-Ninjurin-1 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	(predicted:Human,Mouse,Rat,Chicken,Dog,Cow,Rabbit,Sheep)
Molecular Weight:	Theoretical: 16 kDa.
Purification:	Protein A purified

Applications

Application:	WB,IHC-P,IHC-Fr,ICC/IF,IF,ELISA
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; ICC/IF: 1:100-500; IF: 1:100-500; ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human Ninjurin 1
Antigen Species:	Human
Gene ID:	4814
Uniprot ID:	Q92982
Synonyms:	ninjurin 1
Biology Area:	Epi / Endo-thelial,Neuroregeneration,Neurogenesis,Membrane Proteins

Research Background

Ninjurin family proteins are multi-pass membrane proteins induced by nerve injury in Schwann cells and dorsal root ganglion neurons. Ninjurin proteins act as homophilic cell adhesion molecules that promote axonal growth. Ninjurin proteins also play a role in the formation and function of other tissues. Ninjurin-1 is widely expressed in adult and embryonic tissues, particularly those with epithelial origin. Ninjurin-2 is also widely expressed, with highest levels in adult bone marrow and peripheral blood lymphocytes and embryo liver, thymus and heart. The genes that encode the Ninjurin proteins map to a region known to cause several genetic disorders, including hereditary sensory neuropathy type I and type II (HSN1 and HSN2). However, no link between mutations in the genes encoding Ninjurins and the diseases have been found.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481