

## Anti-Norovirus Capsid protein VP1 Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	(predicted:Norovirus.)
Molecular Weight:	Theoretical: 58 kDa.
Purification:	Protein A purified

### Applications

Application:	ELISA
Recommended	ELISA: 1:5000-10000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	KLH conjugated synthetic peptide: Norovirus Capsid protein VP1
Synonyms:	Norovirus capsid protein;Norwalk virus(strain GI/Human/United States/Norwalk/1968) (Hu/NV/NV/1968/US);p59;Norwalk Viruses;CAPSD_NVN68;NVgp2;Norovirus Capsid protein VP1;58 kd capsid protein [Norwalk virus];VP1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481