


## Anti-PCV2 Cap Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	PCV2
Molecular Weight:	Theoretical: 28 kDa. Actual: 44 kDa.
Purification:	Protein A purified

## Applications

Sample:  
Lane 1: Recombinant PCV2 Cap protein, Trx & His (500ng)  
Verified Activity:  
Primary: Anti-PCV2 Cap protein (TMAB-01346) at 1/1000 dilution  
Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution  
Predicted band size: 28 kDa  
Observed band size: 44 kDa



Application:	WB
Recommended	WB: 1:500-2000

## Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

## Antigen Details

Immunogen:	Recombined PCV2 Cap protein
Synonyms:	cap protein [Porcine circovirus-2]; PCV2 Cap protein; ORF2 capsid protein [ Porcine circovirus 2 ]; CAPSD_PCV2

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481