

Anti-SPNS1 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Rat (predicted:Human,Mouse)
Molecular Weight:	Theoretical: 57 kDa.
Purification:	Protein A purified

Applications

Application:	WB,IHC-P,IHC-Fr,IF,ELISA
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500; ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human SPNS1
Antigen Species:	Human
Gene ID:	83985
Uniprot ID:	Q9H2V7
Synonyms:	SPNS 1;Spinster homolog 3;HSpin1;PP2030;Spinster homolog 1(Drosophila);Protein spinster homolog 1;pp36;Spns3;Spin1;SPINL;LAT;Spinster homolog 2;Spinster;Spinster like;Spns2; Spinster like protein;Spinster-like protein 1;nrs;Spinster homolog 1
Biology Area:	Bcl2 Family

Research Background

Sphingolipid transporter (By similarity). May be involved in necrotic or autophagic cell death.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481