

Anti-TPSB2 Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	(predicted:Human,Mouse,Rat,Pig,Cow,Horse,Rabbit)
Molecular Weight:	Theoretical: 23 kDa.
Purification:	Protein A purified

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: human TPSD1
Antigen Species:	Human
Gene ID:	7177
Uniprot ID:	Q15661
Synonyms:	Tryptase-2;Tryptase II;Tryptase beta-2;TPSB2;TPS2
Biology Area:	Macrophage / Inflamm.,Amino acid metabolism,Cancer,Amino Acids

Research Background

Tryptases comprise a family of trypsin-like serine proteases, the peptidase family S1. Tryptases are enzymatically active only as heparin-stabilized tetramers, and they are resistant to all known endogenous proteinase inhibitors. Several tryptase genes are clustered on chromosome 16p13.3. These genes are characterized by several distinct features. They have a highly conserved 3' UTR and contain tandem repeat sequences at the 5' flank and 3' UTR which are thought to play a role in regulation of the mRNA stability. Although this gene may be an exception, most of the tryptase genes have an intron immediately upstream of the initiator Met codon, which separates the site of transcription initiation from protein coding sequence. This feature is characteristic of tryptases but is unusual in other genes. Tryptases have been implicated as mediators in the pathogenesis of asthma and other allergic and inflammatory disorders. This gene was once considered to be a pseudogene, although it is now believed to be a functional gene that encodes a protein. [provided by RefSeq, Jul 2008]

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481