

Anti-Trophoblast specific protein alpha Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	(predicted:Mouse)
Molecular Weight:	Theoretical: 14 kDa.
Purification:	Protein A purified

Applications

Application:	ELISA,ICC/IF,IF,IHC-Fr,IHC-P,WB
Recommended	WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; ICC/IF: 1:100-500; IF: 1:100-500; ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: mouse Trophoblast specific protein alpha
Antigen Species:	Mouse
Gene ID:	21984
Uniprot ID:	P13438
Synonyms:	Tpbp;TSP_MOUSE;Tpbpa;Tpbpa trophoblast specific protein alpha;Trophoblast specific protein alpha;Trophoblast specific protein;Tb1
Biology Area:	Trophectoderm

Research Background

Trophoblast specific protein alpha is a marker of spongiotrophoblasts.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481