

Anti-Ulex Europaeus Lectin 1 Polyclonal Antibody

Product Details

Ig Type: IgG
Reactivity: (predicted:Ulex europaeus)
Molecular Weight: Theoretical: 27 kDa.
Purification: Protein A purified

Applications

Application: ELISA,ICC/IF,IF,IHC-Fr,IHC-P,WB
Recommended WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; ICC/IF: 1:100-500; IF: 1:100-500; ELISA: 1:5000-10000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: KLH conjugated synthetic peptide: Ulex Europaeus Lectin 1
Synonyms: Anti H(O) lectin I;UEA-1;LEC1_ULEEU;UEA-I;Ulex europaeus agglutinin;Anti H(O) lectin 1;UEA1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481