

Anti-VTG Polyclonal Antibody

Product Details

Ig Type: IgG
Reactivity: (predicted:Medaka)
Purification: Protein A purified

Applications

Application: ELISA
Recommended ELISA: 1:5000-10000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: VTG protein
Synonyms: A0A0F6MTW2_ORYME;vitellogenin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481