

Tween 65

Chemical Properties

CAS No. : 9005-71-4

Formula:

Molecular Weight:

Storage: Pure form: -20°C for 3 years
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Tween 65 is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481