

Sodium carboxymethyl cellulose (MW 250000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sodium carboxymethyl cellulose (CMC-Na) (MW 250000) is the sodium salt of carboxymethyl cellulose. In acidic media, it exhibits adsorption properties and acts as a corrosion inhibitor for low carbon steel. Additionally, sodium carboxymethyl cellulose functions as a thickener, paste, and barrier agent.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481