

## Poloxamer 403 (P123)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at RT

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Poloxamer 403 P123 is a block copolymer of polyethylene oxide and polypropylene oxide, with an average molecular weight of 5750. It exhibits muscle toxicity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481