

Polyglycerol-polyricinoleate

Chemical Properties

CAS No. : 29894-35-7

Formula:

Molecular Weight:

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|-------------------------------------------------------------------------------------------------------------------------|
| Description | Polyglycerol-polyricinoleate (PGPR) is an emulsifier widely used in the food industry to enhance texture and stability. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481