

Dextran T100 (MW 100,000)

Chemical Properties

CAS No. :

Formula: (C₆H₁₀O₅)_n

Molecular Weight:

Storage: Store at RT

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dextran T100 (MW 100,000) is a dehydrated glucose polymer with an average molecular weight of 100,000. This compound exhibits excellent biodegradability and biocompatibility, making it suitable for applications in the food, pharmaceutical, cosmetics, and research sectors.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481