

Dextran T1 (MW 1,000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at RT

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Dextran T1 (MW 1,000) is a polymer of anhydrous glucose with an average molecular weight of 1000. It exhibits excellent biodegradability and biocompatibility, making it suitable for use in the fields of food, pharmaceuticals, cosmetics, and research.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481