

Dextran T3 (MW 3,000)

Chemical Properties

CAS No. :

Formula: (C₆H₁₀O₅)_n

Molecular Weight:

Storage: Store at RT

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dextran T3 (MW 3,000) is a dehydrated glucose polymer with an average molecular weight of 3,000. It exhibits excellent biodegradability and biocompatibility, making it suitable for use in food, pharmaceuticals, cosmetics, and research applications.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481