

DEAE-dextran (MW 40000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight: 40000

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DEAE-dextran (MW 40000) is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.025 mL	0.125 mL	0.250 mL
5 mM	0.005 mL	0.025 mL	0.050 mL
10 mM	0.0025 mL	0.0125 mL	0.025 mL
50 mM	0.0005 mL	0.0025 mL	0.005 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481