

Poly-L-γ-glutamic acid sodium

Chemical Properties

CAS No. : 208106-41-6

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Poly-L-γ-glutamic acid (sodium) is an amino acid polymer.
Targets(IC50)	Amino Acids and Derivatives

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481