

Peptone from soya

Chemical Properties

CAS No. : 91079-46-8

Formula:

Molecular Weight:

Storage: Store at RT

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Peptone from soya (Peptones, soybean) is a type of peptone widely used in culture media. It serves as a nutrient source, providing essential nitrogen, carbon, and other components for microbiology and cell culture. Additionally, Peptone from soya can stimulate or regulate the biosynthesis of cyclic peanoclic acid and exhibits intestinal pool activity in the jejunum of piglets through this mechanism.
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481