

## Iridium(III) chloride hydrate

### Chemical Properties

CAS No. : 14996-61-3

Formula: IrCl<sub>3</sub>·xH<sub>2</sub>O

Molecular Weight:

Storage: Keep away from direct sunlight, Store under nitrogen  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Iridium(III) chloride (hydrate) [Iridium trichloride (hydrate)] serves as a biochemical reagent applicable in the research of biological materials or organic compounds within the field of life sciences.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481