

25-OH Vitamin D3-BSA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

25-OH Vitamin D3-BSA is a hapten-carrier protein conjugate formed by coupling 25-OH Vitamin D3 to bovine serum albumin (BSA). As a hapten, 25-OH Vitamin D3 alone cannot induce an immune response and must be conjugated to a large molecular protein to become a complete antigen with immunogenicity. 25-OH Vitamin D3-BSA may be supplied as a lyophilized powder or as a solution. The solution is prepared in 0.01M PBS (pH 7.4) buffer at a typical concentration of 1mg/mL.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481