

Estradiol-BSA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at -80°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Estradiol-BSA is a hapten-carrier protein conjugate formed by coupling Estradiol to bovine serum albumin (BSA). As a hapten, Estradiol alone cannot induce an immune response and must be conjugated to a large molecular protein to become a complete antigen with immunogenicity. Estradiol-BSA may be supplied as a lyophilized powder or as a solution. The solution is prepared in 0.01M TBS (pH 7.4) buffer at a typical concentration of 1mg/mL.
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481