

## MacConkey Agar

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from moisture,Keep away from direct sunlight

Store at RT

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	MacConkey Agar (MAC) is used for the selective isolation of enteric pathogens. The composition for 1 liter of MacConkey Agar includes: peptone (20.0 g), lactose (10.0 g), bile salts No. 3 (1.5 g), sodium chloride (5.0 g), neutral red (0.03 g), crystal violet (0.001 g), and agar (15.0 g).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481