

Beef Heart Infusion

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from moisture,Keep away from direct sunlight

Store at RT

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Beef Heart Infusion is a non-selective microbial culture medium that supplies essential nutrients for microorganism growth. It is suitable for the cultivation, subculturing, and revitalization of microorganisms.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481