

## Wheat Germ Agglutinin (Fluorescein)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Wheat Germ Agglutinin (WGA) Fluorescein is a plant lectin that specifically binds to biomolecules such as polysaccharides and polypeptides, serving as a probe. It can be used as a biomaterial or organic compound in life sciences research.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481